

Product datasheet for **TL318730**

VCX2 Human shRNA Plasmid Kit (Locus ID 51480)

Product data:

Product Type:	shRNA Plasmids
Product Name:	VCX2 Human shRNA Plasmid Kit (Locus ID 51480)
Locus ID:	51480
Synonyms:	VCX-2r; VCX2R; VCXB
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	VCX2 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 51480). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_016378 , NM_016378.1 , NM_016378.2 , NM_016378.3 , BC096715 , BC098101 , BC098113 , BC098163 , BC104153 , BC104167 , BC104168 , BC134350
UniProt ID:	Q9H322
Summary:	This gene belongs to the VCX/Y gene family, which has multiple members on both X and Y chromosomes that are expressed exclusively in male germ cells. The VCX gene cluster is polymorphic in terms of copy number; different individuals may have a different number of VCX genes. This gene contains two copies of a 30 nt tandem repeat. Deletion of a nearby member of this family was implicated in cognitive disability. [provided by RefSeq, Feb 2015]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).