

Product datasheet for **TL318746**

U2AF1L3 (U2AF1L4) Human shRNA Plasmid Kit (Locus ID 199746)

Product data:

Product Type:	shRNA Plasmids
Product Name:	U2AF1L3 (U2AF1L4) Human shRNA Plasmid Kit (Locus ID 199746)
Locus ID:	199746
Synonyms:	U2AF1-RS3; U2AF1L3; U2AF1L3V1; U2AF1RS3; U2af26
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	U2AF1L4 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 199746). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001040425 , NM_144987 , NR_110171 , NR_110172 , NR_110173 , NM_144987.1 , NM_144987.2 , NM_144987.3 , NM_001040425.1 , NM_001040425.2 , BC021186 , BC021186.1 , BC010865 , BC049208 , BC062474 , BM696851 , NR_163171 , NR_163172 , NR_163173 , NM_001369824 , NR_163169 , NR_163170 , NM_001040425.3
UniProt ID:	Q8WU68
Summary:	RNA-binding protein that function as a pre-mRNA splicing factor. Plays a critical role in both constitutive and enhancer-dependent splicing by mediating protein-protein interactions and protein-RNA interactions required for accurate 3'-splice site selection. Acts by enhancing the binding of U2AF2 to weak pyrimidine tracts. Also participates in the regulation of alternative pre-mRNA splicing. Activates exon 5 skipping of PTPRC during T-cell activation; an event reversed by GFI1. Binds to RNA at the AG dinucleotide at the 3'-splice site (By similarity). Shows a preference for AGC or AGA (By similarity).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).