

## Product datasheet for **TL318884**

### SEC61G Human shRNA Plasmid Kit (Locus ID 23480)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	SEC61G Human shRNA Plasmid Kit (Locus ID 23480)
Locus ID:	23480
Synonyms:	SSS1
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	SEC61G - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 23480). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">BC009480</a> , <a href="#">NM_001012456</a> , <a href="#">NM_014302</a> , <a href="#">NM_014302.1</a> , <a href="#">NM_014302.3</a> , <a href="#">NM_001012456.1</a> , <a href="#">BC009480.2</a> , <a href="#">BC051840</a> , <a href="#">NM_014302.4</a>
UniProt ID:	<a href="#">P60059</a>
Summary:	The Sec61 complex is the central component of the protein translocation apparatus of the endoplasmic reticulum (ER) membrane. Oligomers of the Sec61 complex form a transmembrane channel where proteins are translocated across and integrated into the ER membrane. This complex consists of three membrane proteins- alpha, beta, and gamma. This gene encodes the gamma-subunit protein. Alternatively spliced transcript variants encoding the same protein have been identified. [provided by RefSeq, Jul 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .




[View online »](#)


**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.


For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).

**Product images:**


GFP signal was observed under microscope at 48 hours after transduction of TL318884A virus into HEK293 cells. TL318884A virus was prepared using lenti-shRNA TL318884A and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of TL318884B virus into HEK293 cells. TL318884B virus was prepared using lenti-shRNA TL318884B and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL318884C] virus into HEK293 cells. [TL318884C] virus was prepared using lenti-shRNA [TL318884C] and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL318884D] virus into HEK293 cells. [TL318884D] virus was prepared using lenti-shRNA [TL318884D] and [TR30037] packaging kit.