

## Product datasheet for **TL318923**

### **RNF5 Human shRNA Plasmid Kit (Locus ID 6048)**

#### **Product data:**

Product Type:	shRNA Plasmids
Product Name:	RNF5 Human shRNA Plasmid Kit (Locus ID 6048)
Locus ID:	6048
Synonyms:	RING5; RMA1
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	RNF5 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 6048). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_006913</a> , <a href="#">NM_006913.1</a> , <a href="#">NM_006913.2</a> , <a href="#">NM_006913.3</a> , <a href="#">BC004155</a> , <a href="#">BC020290</a> , <a href="#">BC111392</a> , <a href="#">BC119741</a> , <a href="#">BC119742</a> , <a href="#">BC127651</a> , <a href="#">BC127652</a> , <a href="#">BC148255</a> , <a href="#">BM712189</a> , <a href="#">NM_006913.4</a>
UniProt ID:	<a href="#">Q99942</a>
Summary:	The protein encoded by this gene contains a RING finger, which is a motif known to be involved in protein-protein interactions. This protein is a membrane-bound ubiquitin ligase. It can regulate cell motility by targeting paxillin ubiquitination and altering the distribution and localization of paxillin in cytoplasm and cell focal adhesions. [provided by RefSeq, Jul 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).