

Product datasheet for **TL318927**

RHD Human shRNA Plasmid Kit (Locus ID 6007)

Product data:

Product Type:	shRNA Plasmids
Product Name:	RHD Human shRNA Plasmid Kit (Locus ID 6007)
Locus ID:	6007
Synonyms:	CD240D; DIIIc; RH; Rh4; RH30; RHCED; RhDCw; RHDel; RHDVA(TT); RhII; RhK562-II; RhPI; RHPII; RHXIII
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	RHD - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 6007). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001127691 , NM_001282867 , NM_001282868 , NM_001282869 , NM_001282870 , NM_001282871 , NM_001282872 , NM_016124 , NM_016225 , NM_016124.1 , NM_016124.2 , NM_016124.3 , NM_016124.4 , NM_001127691.1 , NM_001127691.2 , NM_001282867.1 , NM_001282868.1 , NM_001282869.1 , NM_001282872.1 , NM_001282870.1 , NM_001282871.1 , BC139922 , BC160075 , NM_001282871.2 , NM_001127691.3 , NM_001282869.2
UniProt ID:	Q02161
Summary:	The Rh blood group system is the second most clinically significant of the blood groups, second only to ABO. It is also the most polymorphic of the blood groups, with variations due to deletions, gene conversions, and missense mutations. The Rh blood group includes this gene, which encodes the RhD protein, and a second gene that encodes both the RhC and RhE antigens on a single polypeptide. The two genes, and a third unrelated gene, are found in a cluster on chromosome 1. The classification of Rh-positive and Rh-negative individuals is determined by the presence or absence of the highly immunogenic RhD protein on the surface of erythrocytes. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]



[View online »](#)

- shRNA Design:** These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com. If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).
- Performance Guaranteed:** OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.
- For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).