

Product datasheet for **TL319022V**

OBP2A Human shRNA Lentiviral Particle (Locus ID 29991)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	OBP2A Human shRNA Lentiviral Particle (Locus ID 29991)
Locus ID:	29991
Synonyms:	hOBPIIa; LCN13; OBP; OBP2C; OBPIIa
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	OBP2A - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	NM_001293189 , NM_001293193 , NM_014582 , NR_120603 , NM_014582.1 , NM_014582.2 , NM_001293193.1 , NM_001293189.1 , BC069563 , BC069563.1 , BC143803 , NM_001293189.2 , NM_014582.3 , NM_001293193.2
UniProt ID:	Q9NY56
Summary:	This gene encodes a small extracellular protein belonging to the lipocalin superfamily. The protein is thought to transport small, hydrophobic, volatile molecules or odorants through the nasal mucus to olfactory receptors, and may also function as a scavenger of highly concentrated or toxic odors. The protein is expressed as a monomer in the nasal mucus, and can bind diverse types of odorants with a higher affinity for aldehydes and fatty acids. This gene and a highly similar family member are located in a cluster of lipocalin genes on chromosome 9. Alternatively spliced transcript variants have been described, but their biological validity has not been determined. [provided by RefSeq, Jul 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).