

## Product datasheet for **TL319054V**

### NDUFC1 Human shRNA Lentiviral Particle (Locus ID 4717)

#### Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	NDUFC1 Human shRNA Lentiviral Particle (Locus ID 4717)
Locus ID:	4717
Synonyms:	KFYI
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	NDUFC1 - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 <sup>7</sup> TU/ml.
RefSeq:	<a href="#">NM_001184986</a> , <a href="#">NM_001184987</a> , <a href="#">NM_001184988</a> , <a href="#">NM_001184989</a> , <a href="#">NM_001184990</a> , <a href="#">NM_001184991</a> , <a href="#">NM_002494</a> , <a href="#">NM_002494.1</a> , <a href="#">NM_002494.3</a> , <a href="#">NM_001184986.1</a> , <a href="#">NM_001184989.1</a> , <a href="#">NM_001184987.1</a> , <a href="#">NM_001184988.1</a> , <a href="#">NM_001184990.1</a> , <a href="#">NM_001184991.1</a> , <a href="#">BC107682</a> , <a href="#">BC107682.1</a> , <a href="#">BC101798</a> , <a href="#">BC112061</a> , <a href="#">BM768825</a> , <a href="#">NM_001184989.2</a>
UniProt ID:	<a href="#">O43677</a>
Summary:	The encoded protein is a subunit of the NADH:ubiquinone oxidoreductase (complex I), the first enzyme complex in the electron transport chain located in the inner mitochondrial membrane. Alternative splicing results in multiple transcript variants. [provided by RefSeq, May 2010]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).