

Product datasheet for **TL319569**

Galanin (GAL) Human shRNA Plasmid Kit (Locus ID 51083)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Galanin (GAL) Human shRNA Plasmid Kit (Locus ID 51083)
Locus ID:	51083
Synonyms:	ETL8; GAL-GMAP; GALN; GLNN; GMAP
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	GAL - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 51083). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_015973 , NM_015973.1 , NM_015973.2 , NM_015973.3 , NM_015973.4 , BC030241 , BC030241.1 , BM846531
UniProt ID:	P22466
Summary:	This gene encodes a neuroendocrine peptide that is widely expressed in the central and peripheral nervous systems and also the gastrointestinal tract, pancreas, adrenal gland and urogenital tract. The encoded protein is a precursor that is proteolytically processed to generate two mature peptides: galanin and galanin message-associated peptide (GMAP). Galanin has diverse physiological functions including nociception, feeding and energy homeostasis, osmotic regulation and water balance. GMAP has been demonstrated to possess antifungal activity and hypothesized to be part of the innate immune system. [provided by RefSeq, Jul 2015]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .




[View online »](#)


**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.


For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).

Product images:


GFP signal was observed under microscope at 48 hours after transduction of TL319569A virus into HEK293 cells. TL319569A virus was prepared using lenti-shRNA TL319569A and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of TL319569B virus into HEK293 cells. TL319569B virus was prepared using lenti-shRNA TL319569B and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL319569C] virus into HEK293 cells. [TL319569C] virus was prepared using lenti-shRNA [TL319569C] and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL319569D] virus into HEK293 cells. [TL319569D] virus was prepared using lenti-shRNA [TL319569D] and [TR30037] packaging kit.