

Product datasheet for **TL320078**

Bcl2 Binding component 3 (BBC3) Human shRNA Plasmid Kit (Locus ID 27113)

Product data:


Product Type:	shRNA Plasmids
Product Name:	Bcl2 Binding component 3 (BBC3) Human shRNA Plasmid Kit (Locus ID 27113)
Locus ID:	27113
Synonyms:	JFY-1; JFY1; PUMA
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	BBC3 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 27113). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<u>NM_001127240</u> , <u>NM_001127241</u> , <u>NM_001127242</u> , <u>NM_014417</u> , <u>NM_014417.1</u> , <u>NM_014417.2</u> , <u>NM_014417.3</u> , <u>NM_014417.4</u> , <u>NM_001127242.1</u> , <u>NM_001127242.2</u> , <u>NM_001127241.1</u> , <u>NM_001127241.2</u> , <u>NM_001127240.1</u> , <u>NM_001127240.2</u> , <u>BC136481</u> , <u>NM_001127240.3</u> , <u>NM_014417.5</u> , <u>NM_001127242.3</u> , <u>NM_001127241.3</u>
UniProt ID:	<u>Q9BXH1</u>
Summary:	This gene encodes a member of the BCL-2 family of proteins. This family member belongs to the BH3-only pro-apoptotic subclass. The protein cooperates with direct activator proteins to induce mitochondrial outer membrane permeabilization and apoptosis. It can bind to anti-apoptotic Bcl-2 family members to induce mitochondrial dysfunction and caspase activation. Because of its pro-apoptotic role, this gene is a potential drug target for cancer therapy and for tissue injury. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2011]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our <u>custom shRNA service</u> .


[View online »](#)


**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.


For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).

Product images:



GFP signal was observed under microscope at 48 hours after transduction of TL320078A virus into HEK293 cells. TL320078A virus was prepared using lenti-shRNA TL320078A and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of TL320078B virus into HEK293 cells. TL320078B virus was prepared using lenti-shRNA TL320078B and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL320078C] virus into HEK293 cells. [TL320078C] virus was prepared using lenti-shRNA [TL320078C] and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL320078D] virus into HEK293 cells. [TL320078D] virus was prepared using lenti-shRNA [TL320078D] and [TR30037] packaging kit.