

## Product datasheet for **TL320100**

### ALKBH7 Human shRNA Plasmid Kit (Locus ID 84266)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	ALKBH7 Human shRNA Plasmid Kit (Locus ID 84266)
Locus ID:	84266
Synonyms:	ABH7; SPATA11; UNQ6002
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	ALKBH7 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 84266). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_032306</a> , <a href="#">NM_032306.1</a> , <a href="#">NM_032306.2</a> , <a href="#">NM_032306.3</a> , <a href="#">BC004393</a> , <a href="#">BM271861</a> , <a href="#">NM_032306.4</a>
UniProt ID:	<a href="#">Q9BT30</a>
Summary:	May function as protein hydroxylase; can catalyze auto-hydroxylation at Leu-110 (in vitro), but this activity may be due to the absence of the true substrate (PubMed:25122757). Required to induce programmed necrosis in response to DNA damage caused by cytotoxic alkylating agents. Acts by triggering the collapse of mitochondrial membrane potential and loss of mitochondrial function that leads to energy depletion and cell death (PubMed:23666923). ALKBH7-mediated necrosis is probably required to prevent the accumulation of cells with DNA damage (PubMed:23666923). Does not display DNA demethylase activity (PubMed:23666923). Involved in fatty acid metabolism (By similarity). [UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).