

## Product datasheet for **TL320132**

### MEG3 Human shRNA Plasmid Kit (Locus ID 55384)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	MEG3 Human shRNA Plasmid Kit (Locus ID 55384)
Locus ID:	55384
Synonyms:	FP504; GTL2; LINC00023; NCRNA00023; onco-lncRNA-83; prebp1; PRO0518; PRO2160
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	MEG3 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 55384). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NR_002766</a> , <a href="#">NR_003530</a> , <a href="#">NR_003531</a> , <a href="#">NR_033358</a> , <a href="#">NR_033359</a> , <a href="#">NR_033360</a> , <a href="#">NR_046464</a> , <a href="#">NR_046465</a> , <a href="#">NR_046466</a> , <a href="#">NR_046467</a> , <a href="#">NR_046468</a> , <a href="#">NR_046469</a> , <a href="#">NR_046470</a> , <a href="#">NR_046471</a> , <a href="#">NR_046472</a> , <a href="#">NR_046473</a> , <a href="#">BC023543</a> , <a href="#">BC036882</a> , <a href="#">BC051294</a> , <a href="#">BC062783</a> , <a href="#">BC092509</a>
Summary:	This gene is a maternally expressed imprinted gene. Multiple alternatively spliced transcript variants have been transcribed from this gene and all of them are long non-coding RNAs (lncRNAs). This gene is expressed in many normal tissues, but its expression is lost in multiple cancer cell lines of various tissue origins. It inhibits tumor cell proliferation in vitro. It also interacts with the tumor suppressor p53, and regulates p53 target gene expression. Its deletion enhances angiogenesis in vivo. Many experimental evidences demonstrate that this gene is a lncRNA tumor suppressor. [provided by RefSeq, Mar 2012]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .




[View online »](#)


**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).

**Product images:**

GFP signal was observed under microscope at 48 hours after transduction of TL320132B virus into HEK293 cells. TL320132B virus was prepared using lenti-shRNA TL320132B and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL320132D] virus into HEK293 cells. [TL320132D] virus was prepared using lenti-shRNA [TL320132D] and [TR30037] packaging kit.