

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# Product datasheet for TL320259V

## Aryl hydrocarbon Receptor (AHR) Human shRNA Lentiviral Particle (Locus ID 196)

## **Product data:**

Product Type:	shRNA Lentiviral Particles
Product Name:	Aryl hydrocarbon Receptor (AHR) Human shRNA Lentiviral Particle (Locus ID 196)
Locus ID:	196
Synonyms:	bHLHe76; RP85
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	AHR - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10^7 TU/ml.
RefSeq:	<u>NM_001621, NM_001621.1, NM_001621.2, NM_001621.3, NM_001621.4, BC070080, BC070080, BC070080.1, BC021696, BC069390, NM_001621.5</u>
UniProt ID:	<u>P35869</u>
Summary:	The protein encoded by this gene is a ligand-activated helix-loop-helix transcription factor involved in the regulation of biological responses to planar aromatic hydrocarbons. This receptor has been shown to regulate xenobiotic-metabolizing enzymes such as cytochrome P450. Before ligand binding, the encoded protein is sequestered in the cytoplasm; upon ligand binding, this protein moves to the nucleus and stimulates transcription of target genes. [provided by RefSeq, Sep 2015]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <u>techsupport@origene.com</u> . If you need a special design or shRNA sequence, please utilize our <u>custom shRNA service</u> .




#### CRIGENE Aryl hydrocarbon Receptor (AHR) Human shRNA Lentiviral Particle (Locus ID 196) – TL320259V


Performance Guaranteed: OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).

### **Product images:**




GFP signal was observed under microscope at 48 hours after transduction of TL320259A virus into HEK293 cells. TL320259A virus was prepared using lenti-shRNA TL320259A and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of TL320259A virus into HEK293 cells. TL320259A virus was prepared using lenti-shRNA TL320259A and [TR30037] packaging kit.

GFP signal was observed under microscope at 48 hours after transduction of TL320259B virus into HEK293 cells. TL320259B virus was prepared using lenti-shRNA TL320259B and [TR30037] packaging kit.


GFP signal was observed under microscope at 48 hours after transduction of TL320259B virus into HEK293 cells. TL320259B virus was prepared using lenti-shRNA TL320259B and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL320259C] virus into HEK293 cells. [TL320259C] virus was prepared using lenti-shRNA [TL320259C] and [TR30037] packaging kit.

GFP signal was observed under microscope at 48 hours after transduction of [TL320259C] virus into HEK293 cells. [TL320259C] virus was prepared using lenti-shRNA [TL320259C] and [TR30037] packaging kit.

GFP signal was observed under microscope at 48 hours after transduction of [TL320259D] virus into HEK293 cells. [TL320259D] virus was prepared using lenti-shRNA [TL320259D] and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL320259D] virus into HEK293 cells. [TL320259D] virus was prepared using lenti-shRNA [TL320259D] and [TR30037] packaging kit.