

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TL320317

CKII alpha (CSNK2A1) Human shRNA Plasmid Kit (Locus ID 1457)

Product data:

Product Type:	shRNA Plasmids
Product Name:	CKII alpha (CSNK2A1) Human shRNA Plasmid Kit (Locus ID 1457)
Locus ID:	1457
Synonyms:	CK2A1; Cka1; Cka2; CKII; OCNDS
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	CSNK2A1 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 1457). 5μg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM 001895, NM 177559, NM 177560, NM 177560.1, NM 177560.2, NM 177559.1, NM 177559.2, NM 001895.1, NM 001895.2, NM 001895.3, BC053532, BC053532.1, BC050036, BC002615, BC011668, BC071167, NM 001362771, NM 001362770, NM 001895.4, NM 177559.3
UniProt ID:	<u>P68400</u>
Summary:	Casein kinase II is a serine/threonine protein kinase that phosphorylates acidic proteins such as casein. It is involved in various cellular processes, including cell cycle control, apoptosis, and circadian rhythm. The kinase exists as a tetramer and is composed of an alpha, an alpha- prime, and two beta subunits. The alpha subunits contain the catalytic activity while the beta subunits undergo autophosphorylation. The protein encoded by this gene represents the alpha subunit. Multiple transcript variants encoding different protein isoforms have been found for this gene. [provided by RefSeq, Apr 2018]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <u>techsupport@origene.com</u> . If you need a special design or shRNA sequence, please utilize our <u>custom shRNA service</u> .



View online »


This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

GRIGENE CKII alpha (CSNK2A1) Human shRNA Plasmid Kit (Locus ID 1457) – TL320317


Performance Guaranteed: OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).

Product images:



GFP signal was observed under microscope at 48 hours after transduction of TL320317A virus into HEK293 cells. TL320317A virus was prepared using lenti-shRNA TL320317A and [TR30037] packaging kit.

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US 

GFP signal was observed under microscope at 48 hours after transduction of TL320317B virus into HEK293 cells. TL320317B virus was prepared using lenti-shRNA TL320317B and [TR30037] packaging kit.

GFP signal was observed under microscope at 48 hours after transduction of [TL320317C] virus into HEK293 cells. [TL320317C] virus was prepared using lenti-shRNA [TL320317C] and [TR30037] packaging kit.

GFP signal was observed under microscope at 48 hours after transduction of [TL320317D] virus into HEK293 cells. [TL320317D] virus was prepared using lenti-shRNA [TL320317D] and [TR30037] packaging kit.

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US