

## Product datasheet for **TL320359V**

### FOXO3 Human shRNA Lentiviral Particle (Locus ID 2309)

#### Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	FOXO3 Human shRNA Lentiviral Particle (Locus ID 2309)
Locus ID:	2309
Synonyms:	AF6q21; FKHRL1; FKHRL1P2; FOXO2; FOXO3A
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	FOXO3 - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 <sup>7</sup> TU/ml.
RefSeq:	<a href="#">NM_001455</a> , <a href="#">NM_201559</a> , <a href="#">NM_201559.1</a> , <a href="#">NM_201559.2</a> , <a href="#">NM_001455.1</a> , <a href="#">NM_001455.2</a> , <a href="#">NM_001455.3</a> , <a href="#">BC021224</a> , <a href="#">BC021224.2</a> , <a href="#">BC068552</a> , <a href="#">BC020227</a> , <a href="#">BC045800</a> , <a href="#">NM_201559.3</a>
UniProt ID:	<a href="#">O43524</a>
Summary:	This gene belongs to the forkhead family of transcription factors which are characterized by a distinct forkhead domain. This gene likely functions as a trigger for apoptosis through expression of genes necessary for cell death. Translocation of this gene with the MLL gene is associated with secondary acute leukemia. Alternatively spliced transcript variants encoding the same protein have been observed. [provided by RefSeq, Jul 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).