

## Product datasheet for **TL320430**

### NEK3 Human shRNA Plasmid Kit (Locus ID 4752)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	NEK3 Human shRNA Plasmid Kit (Locus ID 4752)
Locus ID:	4752
Synonyms:	HSPK36
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	NEK3 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 4752). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001146099</a> , <a href="#">NM_002498</a> , <a href="#">NM_152720</a> , <a href="#">NR_027415</a> , <a href="#">NM_152720.1</a> , <a href="#">NM_152720.2</a> , <a href="#">NM_002498.1</a> , <a href="#">NM_002498.2</a> , <a href="#">NM_001146099.1</a> , <a href="#">BC019916</a> , <a href="#">BC019916.1</a> , <a href="#">BM675414</a> , <a href="#">BM984985</a> , <a href="#">NM_002498.3</a>
UniProt ID:	<a href="#">P51956</a>
Summary:	This gene encodes a member of the NimA (never in mitosis A) family of serine/threonine protein kinases. The encoded protein differs from other NimA family members in that it is not cell cycle regulated and is found primarily in the cytoplasm. The kinase is activated by prolactin stimulation, leading to phosphorylation of VAV2 guanine nucleotide exchange factor, paxillin, and activation of the RAC1 GTPase. Two functional alleles for this gene have been identified in humans. The reference genome assembly (GRCh38) represents a functional allele that is associated with the inclusion of an additional coding exon in protein-coding transcripts, compared to an alternate functional allele that lacks the exon. [provided by RefSeq, Sep 2019]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).