

## Product datasheet for **TL320440V**

### DDR2 Human shRNA Lentiviral Particle (Locus ID 4921)

#### Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	DDR2 Human shRNA Lentiviral Particle (Locus ID 4921)
Locus ID:	4921
Synonyms:	MIG20a; NTRKR3; TKT; TYRO10; WRCN
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	DDR2 - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 <sup>7</sup> TU/ml.
RefSeq:	<a href="#">NM_001014796</a> , <a href="#">NM_006182</a> , <a href="#">NM_001354982</a> , <a href="#">NM_001354983</a> , <a href="#">NM_001014796.1</a> , <a href="#">NM_006182.1</a> , <a href="#">NM_006182.2</a> , <a href="#">BC052998</a> , <a href="#">BC052998.2</a> , <a href="#">NM_001014796.3</a> , <a href="#">NM_006182.4</a>
UniProt ID:	<a href="#">Q16832</a>
Summary:	This gene encodes a member of the discoidin domain receptor subclass of the receptor tyrosine kinase (RTKs) protein family. RTKs play a key role in the communication of cells with their microenvironment. The encoded protein is a collagen-induced receptor that activates signal transduction pathways involved in cell adhesion, proliferation, and extracellular matrix remodeling. This protein is expressed in numerous cell types and may also be involved in wound repair and regulate tumor growth and invasiveness. Mutations in this gene are the cause of short limb-hand type spondylometaphyseal dysplasia. [provided by RefSeq, Aug 2017]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .




[View online »](#)


**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.


For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).

**Product images:**


GFP signal was observed under microscope at 48 hours after transduction of TL320440A virus into HEK293 cells. TL320440A virus was prepared using lenti-shRNA TL320440A and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of TL320440B virus into HEK293 cells. TL320440B virus was prepared using lenti-shRNA TL320440B and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL320440C] virus into HEK293 cells. [TL320440C] virus was prepared using lenti-shRNA [TL320440C] and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL320440D] virus into HEK293 cells. [TL320440D] virus was prepared using lenti-shRNA [TL320440D] and [TR30037] packaging kit.