

Product datasheet for **TL320546V**

TIE2 (TEK) Human shRNA Lentiviral Particle (Locus ID 7010)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	TIE2 (TEK) Human shRNA Lentiviral Particle (Locus ID 7010)
Locus ID:	7010
Synonyms:	CD202B; GLC3E; TIE-2; TIE2; VMCM; VMCM1
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	TEK - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	NM_000459 , NM_001290077 , NM_001290078 , NM_000459.1 , NM_000459.2 , NM_000459.3 , NM_000459.4 , NM_001290078.1 , NM_001290077.1 , BC035514 , BC035514.1 , NM_000459.5
UniProt ID:	Q02763
Summary:	This gene encodes a receptor that belongs to the protein tyrosine kinase Tie2 family. The encoded protein possesses a unique extracellular region that contains two immunoglobulin-like domains, three epidermal growth factor (EGF)-like domains and three fibronectin type III repeats. The ligand angiopoietin-1 binds to this receptor and mediates a signaling pathway that functions in embryonic vascular development. Mutations in this gene are associated with inherited venous malformations of the skin and mucous membranes. Alternative splicing results in multiple transcript variants. Additional alternatively spliced transcript variants of this gene have been described, but their full-length nature is not known. [provided by RefSeq, Feb 2014]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .




[View online »](#)


**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.


For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).

Product images:


GFP signal was observed under microscope at 48 hours after transduction of TL320546A virus into HEK293 cells. TL320546A virus was prepared using lenti-shRNA TL320546A and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of TL320546B virus into HEK293 cells. TL320546B virus was prepared using lenti-shRNA TL320546B and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL320546C] virus into HEK293 cells. [TL320546C] virus was prepared using lenti-shRNA [TL320546C] and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL320546D] virus into HEK293 cells. [TL320546D] virus was prepared using lenti-shRNA [TL320546D] and [TR30037] packaging kit.