

Product datasheet for **TL320557**

TNFAIP3 Human shRNA Plasmid Kit (Locus ID 7128)

Product data:

Product Type:	shRNA Plasmids
Product Name:	TNFAIP3 Human shRNA Plasmid Kit (Locus ID 7128)
Locus ID:	7128
Synonyms:	A20; AISBL; OTUD7C; TNFA1P2
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	TNFAIP3 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 7128). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001270507 , NM_001270508 , NM_006290 , NM_006290.1 , NM_006290.2 , NM_006290.3 , NM_001270507.1 , NM_001270508.1 , BC041790 , BC064689 , BC098379 , BC113871 , BC114480 , NM_001270508.2 , NM_001270507.2 , NM_006290.4
UniProt ID:	P21580
Summary:	This gene was identified as a gene whose expression is rapidly induced by the tumor necrosis factor (TNF). The protein encoded by this gene is a zinc finger protein and ubiquitin-editing enzyme, and has been shown to inhibit NF-kappa B activation as well as TNF-mediated apoptosis. The encoded protein, which has both ubiquitin ligase and deubiquitinase activities, is involved in the cytokine-mediated immune and inflammatory responses. Several transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Jul 2012]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .




[View online »](#)


**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.


For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).

Product images:


GFP signal was observed under microscope at 48 hours after transduction of TL320557A virus into HEK293 cells. TL320557A virus was prepared using lenti-shRNA TL320557A and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of TL320557B virus into HEK293 cells. TL320557B virus was prepared using lenti-shRNA TL320557B and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL320557C] virus into HEK293 cells. [TL320557C] virus was prepared using lenti-shRNA [TL320557C] and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL320557D] virus into HEK293 cells. [TL320557D] virus was prepared using lenti-shRNA [TL320557D] and [TR30037] packaging kit.