

## Product datasheet for **TL320609**

### MAPKAP Kinase 2 (MAPKAPK2) Human shRNA Plasmid Kit (Locus ID 9261)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	MAPKAP Kinase 2 (MAPKAPK2) Human shRNA Plasmid Kit (Locus ID 9261)
Locus ID:	9261
Synonyms:	MAPKAP-K2; MK-2; MK2
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	MAPKAPK2 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 9261). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_004759</a> , <a href="#">NM_032960</a> , <a href="#">NM_032960.1</a> , <a href="#">NM_032960.2</a> , <a href="#">NM_032960.3</a> , <a href="#">NM_004759.1</a> , <a href="#">NM_004759.2</a> , <a href="#">NM_004759.3</a> , <a href="#">NM_004759.4</a> , <a href="#">BC052584</a> , <a href="#">BC036060</a> , <a href="#">NM_032960.4</a> , <a href="#">NM_004759.5</a>
UniProt ID:	<a href="#">P49137</a>
Summary:	This gene encodes a member of the Ser/Thr protein kinase family. This kinase is regulated through direct phosphorylation by p38 MAP kinase. In conjunction with p38 MAP kinase, this kinase is known to be involved in many cellular processes including stress and inflammatory responses, nuclear export, gene expression regulation and cell proliferation. Heat shock protein HSP27 was shown to be one of the substrates of this kinase in vivo. Two transcript variants encoding two different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).