

## Product datasheet for **TL320641**

### STK25 Human shRNA Plasmid Kit (Locus ID 10494)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	STK25 Human shRNA Plasmid Kit (Locus ID 10494)
Locus ID:	10494
Synonyms:	SOK1; YSK1
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	STK25 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 10494). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001271977</a> , <a href="#">NM_001271978</a> , <a href="#">NM_001271979</a> , <a href="#">NM_001271980</a> , <a href="#">NM_001282305</a> , <a href="#">NM_001282306</a> , <a href="#">NM_001282307</a> , <a href="#">NM_001282308</a> , <a href="#">NM_006374</a> , <a href="#">NR_073530</a> , <a href="#">NR_073531</a> , <a href="#">NR_073532</a> , <a href="#">NR_073533</a> , <a href="#">NM_006374.1</a> , <a href="#">NM_006374.2</a> , <a href="#">NM_006374.3</a> , <a href="#">NM_006374.4</a> , <a href="#">NM_001271979.1</a> , <a href="#">NM_001271980.1</a> , <a href="#">NM_001271977.1</a> , <a href="#">NM_001271978.1</a> , <a href="#">NM_001282305.1</a> , <a href="#">NM_001282307.1</a> , <a href="#">NM_001282308.1</a> , <a href="#">NM_001282306.1</a> , <a href="#">BC007852</a> , <a href="#">BC007852.2</a> , <a href="#">BC015793</a> , <a href="#">BC091505</a> , <a href="#">NM_001271979.2</a> , <a href="#">NM_001271977.2</a> , <a href="#">NM_001271978.2</a> , <a href="#">NM_001282306.2</a> , <a href="#">NM_006374.5</a> , <a href="#">NM_001271980.2</a> , <a href="#">NM_001282308.2</a>
UniProt ID:	<a href="#">O00506</a>
Summary:	This gene encodes a member of the germinal centre kinase III (GCK III) subfamily of the sterile 20 superfamily of kinases. The encoded enzyme plays a role in serine-threonine liver kinase B1 (LKB1) signaling pathway to regulate neuronal polarization and morphology of the Golgi apparatus. The protein is translocated from the Golgi apparatus to the nucleus in response to chemical anoxia and plays a role in regulation of cell death. A pseudogene associated with this gene is located on chromosome 18. Multiple alternatively spliced transcript variants have been observed for this gene. [provided by RefSeq, Dec 2012]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).