

Product datasheet for **TL320934**

AKR1B15 Human shRNA Plasmid Kit (Locus ID 441282)

Product data:

Product Type:	shRNA Plasmids
Product Name:	AKR1B15 Human shRNA Plasmid Kit (Locus ID 441282)
Locus ID:	441282
Synonyms:	AK1R1B7; AKR1B10L; AKR1R1B7
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	AKR1B15 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 441282). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001080538 , NM_001080538.1 , NM_001080538.2 , BC072676 , NM_001367820 , NM_001367821 , NR_160301 , NR_160302 , NM_001367822 , NM_001080538.3
UniProt ID:	C9JRZ8
Summary:	Isoform 1: Catalyzes the NADPH-dependent reduction of a variety of carbonyl substrates, like aromatic aldehydes, alkenals, ketones and alpha-dicarbonyl compounds (PubMed:26222439, PubMed:21276782). In addition, catalyzes the reduction of androgens and estrogens with high positional selectivity (shows 17-beta-hydroxysteroid dehydrogenase activity) as well as 3-keto-acyl-CoAs (PubMed:25577493). Displays strong enzymatic activity toward all-trans-retinal and 9-cis-retinal (PubMed:26222439). May play a physiological role in retinoid metabolism (PubMed:26222439).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).