

Product datasheet for **TL320972**

LYPD8 Human shRNA Plasmid Kit (Locus ID 646627)

Product data:

Product Type:	shRNA Plasmids
Product Name:	LYPD8 Human shRNA Plasmid Kit (Locus ID 646627)
Locus ID:	646627
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	LYPD8 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 646627). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001085474 , NM_001291283 , NM_001085474.1 , NM_001291283.1 , BC033943
UniProt ID:	Q6UX82
Summary:	Secreted protein specifically required to prevent invasion of Gram-negative bacteria in the inner mucus layer of the colon epithelium, a portion of the large intestine which is free of commensal microbiota. Prevents invasion of flagellated microbiota by binding to the flagellum of bacteria, such as P.mirabilis, thereby inhibiting bacterial motility in the intestinal lumen. Segregation of intestinal bacteria and epithelial cells in the colon is required to preserve intestinal homeostasis.[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).