

Product datasheet for **TL321062**

RSPH10B2 Human shRNA Plasmid Kit (Locus ID 728194)

Product data:

Product Type:	shRNA Plasmids
Product Name:	RSPH10B2 Human shRNA Plasmid Kit (Locus ID 728194)
Locus ID:	728194
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	RSPH10B2 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 728194). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC044242 , NM_001099697 , NM_001099697.1 , BC044242.1
UniProt ID:	B2RC85
Summary:	This gene encodes a protein component of the radial spoke head in flagella and motile cilia. Eukaryotic flagella and motile cilia share a common 9 + 2 structure, in which nine peripheral microtubule doublets (MTDs) surround a central-pair of microtubules (CP), with radial spokes connecting the MTDs to the CP. The radial spoke is a multi-protein complex that works as a mechanochemical transducer between the CP and the MTDs. The radial spoke contributes to the regulation of the activity of dynein motors, and thus to flagellar motility. PMID: 22754630 provides a good review of radial spokes. [provided by RefSeq, Jul 2017]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .


[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).