

Product datasheet for **TL32112V**

MARCHF11 Human shRNA Lentiviral Particle (Locus ID 441061)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	MARCHF11 Human shRNA Lentiviral Particle (Locus ID 441061)
Locus ID:	441061
Synonyms:	MARCH-XI; MARCH11; RNF226
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	MARCH11 - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	NM_001102562 , NM_001102562.1 , BC150513 , BC150514 , NM_001102562.3
UniProt ID:	A6NNE9
Summary:	MARCH11 is a member of the MARCH family of membrane-bound E3 ubiquitin ligases (EC 6.3.2.19). These enzymes add ubiquitin (see MIM 191339) to target lysines in substrate proteins, thereby signaling their intracellular transport. March11 appears to have a role in ubiquitin-mediated protein sorting in the trans-Golgi network (TGN)-multivesicular body (MVB) transport pathway (Morokuma et al., 2007 [PubMed 17604280]).[supplied by OMIM, Apr 2010]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).