

## Product datasheet for **TL321487**

### CECR2 Human shRNA Plasmid Kit (Locus ID 27443)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	CECR2 Human shRNA Plasmid Kit (Locus ID 27443)
Locus ID:	27443
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	CECR2 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 27443). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001290046</a> , <a href="#">NM_001290047</a> , <a href="#">NM_031413</a> , <a href="#">NM_031413.1</a> , <a href="#">NM_031413.2</a> , <a href="#">NM_031413.3</a> , <a href="#">NM_031413.4</a> , <a href="#">NM_001290046.1</a> , <a href="#">NM_001290047.1</a> , <a href="#">BC166664</a> , <a href="#">NM_001290047.2</a>
UniProt ID:	<a href="#">Q9BXF3</a>
Summary:	This gene encodes a bromodomain-containing protein that is involved in chromatin remodeling, and may additionally play a role in DNA damage response. The encoded protein functions as part of an ATP-dependent complex that is involved in neurulation. This gene is a candidate gene for Cat Eye Syndrome. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2014]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .




[View online »](#)


**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.


For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).

**Product images:**


GFP signal was observed under microscope at 48 hours after transduction of TL321487A virus into HEK293 cells. TL321487A virus was prepared using lenti-shRNA TL321487A and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of TL321487B virus into HEK293 cells. TL321487B virus was prepared using lenti-shRNA TL321487B and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL321487C] virus into HEK293 cells. [TL321487C] virus was prepared using lenti-shRNA [TL321487C] and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL321487D] virus into HEK293 cells. [TL321487D] virus was prepared using lenti-shRNA [TL321487D] and [TR30037] packaging kit.