

Product datasheet for **TL500009**

Acads Mouse shRNA Plasmid (Locus ID 11409)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Acads Mouse shRNA Plasmid (Locus ID 11409)
Locus ID:	11409
Synonyms:	AI196007; Bcd-1; Bcd1; Hdlq; Hdlq8; SC; SCAD
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Acads - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 11409). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC016259 , NM_007383 , NM_007383.1 , NM_007383.2 , NM_007383.3 , BC055063
UniProt ID:	Q07417
Summary:	This gene encodes a homotetrameric mitochondrial flavoprotein and is a member of the acyl-CoA dehydrogenase family. Members of this family catalyze the first step of fatty acid beta-oxidation, forming a C2-C3 trans-double bond in a FAD-dependent reaction. As beta-oxidation cycles through its four steps, each member of the Acyl-CoA dehydrogenase family works at an optimum fatty acid chain-length. This enzyme has its optimum at C(four)-CoA. In mice, deficiency of this gene has been linked to cold sensitivity and increased high-density lipoprotein levels. [provided by RefSeq, Nov 2012]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).