

## Product datasheet for **TL500041**

### Adamts1 Mouse shRNA Plasmid (Locus ID 11504)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Adamts1 Mouse shRNA Plasmid (Locus ID 11504)
Locus ID:	11504
Synonyms:	ADAM-; ADAM-TS1; ADAMT; ADAMTS; ADAMTS-1; C3-C5; METH-1; METH1
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Adamts1 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 11504). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">BC040382</a> , <a href="#">BC050834</a> , <a href="#">NM_009621</a> , <a href="#">NM_009621.1</a> , <a href="#">NM_009621.2</a> , <a href="#">NM_009621.3</a> , <a href="#">NM_009621.4</a> , <a href="#">NM_009621.5</a> , <a href="#">BC003269</a>
UniProt ID:	<a href="#">P97857</a>
Summary:	This gene encodes a member of the ADAMTS (a disintegrin and metalloproteinase with thrombospondin motif) family and preproprotein that is proteolytically processed to generate a mature protein product. This secreted protein product plays an important role in ovulation, likely through its cleavage of the extracellular matrix component versican. The encoded protein may enhance tumorigenesis in a mouse model of breast cancer. Homozygous knockout mice for this gene exhibit enhanced perinatal lethality, impaired growth and adipose tissue development, and impaired ovulation in females. [provided by RefSeq, Oct 2015]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).