

## Product datasheet for **TL500257**

### Camk2g Mouse shRNA Plasmid (Locus ID 12325)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Camk2g Mouse shRNA Plasmid (Locus ID 12325)
Locus ID:	12325
Synonyms:	Camkg
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Camk2g - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 12325). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">BC019162</a> , <a href="#">BC025597</a> , <a href="#">NM_001039138</a> , <a href="#">NM_001039139</a> , <a href="#">NM_178597</a> , <a href="#">NM_001039139.1</a> , <a href="#">NM_001039139.2</a> , <a href="#">NM_001039138.1</a> , <a href="#">NM_001039138.2</a> , <a href="#">NM_178597.1</a> , <a href="#">NM_178597.2</a> , <a href="#">NM_178597.3</a> , <a href="#">NM_178597.4</a> , <a href="#">NM_178597.5</a> , <a href="#">NM_001360186</a>
UniProt ID:	<a href="#">Q923T9</a>
Summary:	Calcium/calmodulin-dependent protein kinase that functions autonomously after Ca(2+)/calmodulin-binding and autophosphorylation, and is involved in sarcoplasmic reticulum Ca(2+) transport in skeletal muscle and may function in dendritic spine and synapse formation and neuronal plasticity. In slow-twitch muscles, is involved in regulation of sarcoplasmic reticulum (SR) Ca(2+) transport and in fast-twitch muscle participates in the control of Ca(2+) release from the SR through phosphorylation of the ryanodine receptor-coupling factor triadin. In neurons, may participate in the promotion of dendritic spine and synapse formation and maintenance of synaptic plasticity which enables long-term potentiation (LTP) and hippocampus-dependent learning (By similarity).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).