

## Product datasheet for **TL500494V**

### Dach1 Mouse shRNA Lentiviral Particle (Locus ID 13134)

#### Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Dach1 Mouse shRNA Lentiviral Particle (Locus ID 13134)
Locus ID:	13134
Synonyms:	Dac; Dach
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	Dach1 - Mouse shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 <sup>7</sup> TU/ml.
RefSeq:	<a href="#">NM_001038610</a> , <a href="#">NM_007826</a> , <a href="#">NM_007826.1</a> , <a href="#">NM_007826.2</a> , <a href="#">NM_007826.3</a> , <a href="#">NM_001038610.1</a> , <a href="#">NM_001038610.2</a> , <a href="#">BC141130</a> , <a href="#">BC052005</a> , <a href="#">BC059934</a> , <a href="#">BC078644</a> , <a href="#">BC145440</a>
UniProt ID:	<a href="#">Q9QYB2</a>
Summary:	Transcription factor that is involved in regulation of organogenesis. Seems to be a regulator of SIX1, SIX6 and probably SIX5. Corepression of precursor cell proliferation in myoblasts by SIX1 is switched to coactivation through recruitment of EYA3 to the SIX1-DACH1 complex. Transcriptional activation seems also to involve association of CREBBP. Seems to act as a corepressor of SIX6 in regulating proliferation by directly repressing cyclin-dependent kinase inhibitors, including the p27Kip1 promoter. Inhibits TGF-beta signaling through interaction with SMAD4 and NCOR1 (By similarity). Binds to chromatin DNA via its DACHbox-N domain. [UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).