

## Product datasheet for **TL500669**

### Srsf10 Mouse shRNA Plasmid (Locus ID 14105)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Srsf10 Mouse shRNA Plasmid (Locus ID 14105)
Locus ID:	14105
Synonyms:	Fusip1; FUSIP2; Nssr; NSSR1; NSSR2; Sfrs13a; SRrp40; Srsf13a; TASR; TASR1; TASR2
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Srsf10 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 14105). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">BC037591</a> , <a href="#">BC043060</a> , <a href="#">BC083082</a> , <a href="#">NM_001080387</a> , <a href="#">NM_001284195</a> , <a href="#">NM_001284196</a> , <a href="#">NM_010178</a> , <a href="#">NM_001080387.1</a> , <a href="#">NM_001080387.2</a> , <a href="#">NM_010178.1</a> , <a href="#">NM_010178.2</a> , <a href="#">NM_010178.3</a> , <a href="#">NM_001284196.1</a> , <a href="#">NM_001284195.1</a> , <a href="#">BC037591.1</a> , <a href="#">BC006670</a>
UniProt ID:	<a href="#">Q9R0U0</a>
Summary:	Splicing factor that in its dephosphorylated form acts as a general repressor of pre-mRNA splicing. Seems to interfere with the U1 snRNP 5'-splice recognition of SNRNP70. Required for splicing repression in M-phase cells and after heat shock. Also acts as a splicing factor that specifically promotes exon skipping during alternative splicing. Interaction with YTHDC1, a RNA-binding protein that recognizes and binds N6-methyladenosine (m6A)-containing RNAs, prevents SRSF10 from binding to its mRNA-binding sites close to m6A-containing regions, leading to inhibit exon skipping during alternative splicing (By similarity). May be involved in regulation of alternative splicing in neurons (PubMed:10583508).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).