

Product datasheet for **TL500936**

Hcfc1 Mouse shRNA Plasmid (Locus ID 15161)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Hcfc1 Mouse shRNA Plasmid (Locus ID 15161)
Locus ID:	15161
Synonyms:	AW060991; HCF; HCF1
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Hcfc1 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 15161). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC053742 , NM_008224 , NM_008224.1 , NM_008224.2 , NM_008224.3 , NM_008224.4 , BC030898
UniProt ID:	Q61191
Summary:	This gene encodes a transcription cofactor that regulates the progression of cell cycle and maintain the ability of embryonic stem cells to self-renew. The encoded protein is a large precursor that undergoes site-specific proteolytic cleavage to yield N- and C-terminal chains that form a non-covalent heterodimer. The encoded protein has been implicated in the regulation of expression of immediate early genes after herpes simplex virus infection and glucose-stimulated secretion of insulin by pancreatic beta cells. [provided by RefSeq, Aug 2015]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).