

Product datasheet for **TL501064V**

Cd74 Mouse shRNA Lentiviral Particle (Locus ID 16149)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Cd74 Mouse shRNA Lentiviral Particle (Locus ID 16149)
Locus ID:	16149
Synonyms:	CLIP; DHLAG; HLADG; Ia-GAMMA; li
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	Cd74 - Mouse shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, $>10^7$ TU/ml.
RefSeq:	BC003476 , BC061489 , BC096435 , NM_001042605 , NM_010545 , NM_010545.1 , NM_010545.2 , NM_010545.3 , NM_001042605.1
UniProt ID:	P04441
Summary:	Plays a critical role in MHC class II antigen processing by stabilizing peptide-free class II alpha/beta heterodimers in a complex soon after their synthesis and directing transport of the complex from the endoplasmic reticulum to compartments where peptide loading of class II takes place. Enhance also the stimulation of T-cell responses through interaction with CD44.[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .

[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).