

Product datasheet for **TL501124**

Itih1 Mouse shRNA Plasmid (Locus ID 16424)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Itih1 Mouse shRNA Plasmid (Locus ID 16424)
Locus ID:	16424
Synonyms:	In; Intin1; ITI-HC1; Itih; Itih-1
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Itih1 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 16424). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC028814 , NM_008406 , NM_008406.1 , NM_008406.2 , NM_008406.3 , BC028814.1 , BC013465
UniProt ID:	Q61702
Summary:	This gene encodes a heavy chain of inter-alpha trypsin inhibitor (IaI) family of plasma serine protease inhibitors. IaI proteins are protein-glycosaminoglycan-protein complexes comprised of two heavy chains and a light chain. The encoded protein covalently associates with the light chain via a chondroitin sulfate moiety. Intravenous administration of the encoded protein improved survival of mice after infection with Escherichia coli. This gene is located adjacent to two other IaI heavy chain genes. Alternative splicing results in multiple transcript variants encoding different isoforms that may undergo similar proteolytic processing to generate mature protein. [provided by RefSeq, Oct 2015]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).