

Product datasheet for **TL501195**

Klk1b26 Mouse shRNA Plasmid (Locus ID 16618)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Klk1b26 Mouse shRNA Plasmid (Locus ID 16618)
Locus ID:	16618
Synonyms:	Egfbp; Egfbp2; Klk2; Klk26; mGK-26; mGK26; PRECE; PRECE-2
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Klk1b26 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 16618). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC024688 , NM_010644 , NM_010644.1 , NM_010644.2 , NM_010644.3 , BC138704
UniProt ID:	P36369
Summary:	This gene encodes a member of the kallikrein subfamily of the peptidase S1 family of serine proteases. The kallikrein genes are present in a large gene cluster on chromosome 7. The encoded preproprotein is proteolytically processed to generate a mature protein product. Expression of this gene in the submandibular gland is sexually dimorphic likely due to both transcriptional and post-transcriptional regulation. This gene is thought to be distinct from the Egfbp2 gene (Gene ID: 13647), with which it shares 98% identity (PMIDs: 1959648, 9685728), however, it is not clear if both genes are present in all strains of mice. [provided by RefSeq, Sep 2015]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .


[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).