

Product datasheet for TL501209V

Kras Mouse shRNA Lentiviral Particle (Locus ID 16653)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Kras Mouse shRNA Lentiviral Particle (Locus ID 16653)
Locus ID:	16653
Synonyms:	Al929937; c-K-ras; c-Ki-ras; K-Ras; K-ras; Ki-ras; Kras-2; Kras2; p21B; ras
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	Kras - Mouse shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10^7 TU/ml.
RefSeq:	<u>BC004642, BC010202, BC119495, NM 021284, NM 021284.2, NM 021284.3, NM 021284.4, NM 021284.5, NM 021284.6, BC014700, BC031623, BC141051</u>
UniProt ID:	<u>P32883</u>
Summary:	Ras proteins bind GDP/GTP and possess intrinsic GTPase activity (By similarity). Plays an important role in the regulation of cell proliferation (PubMed:6474169, PubMed:1352876). Plays a role in promoting oncogenic events by inducing transcriptional silencing of tumor suppressor genes (TSGs) in colorectal cancer (CRC) cells in a ZNF304-dependent manner (By similarity).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <u>techsupport@origene.com</u> . If you need a special design or shRNA sequence, please utilize our custom shRNA service.

View online »

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

GRIGENE Kras Mouse shRNA Lentiviral Particle (Locus ID 16653) – TL501209V

Performance Guaranteed: OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US