

Product datasheet for **TL501496V**

Notch2 Mouse shRNA Lentiviral Particle (Locus ID 18129)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Notch2 Mouse shRNA Lentiviral Particle (Locus ID 18129)
Locus ID:	18129
Synonyms:	AI853703; N2
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	Notch2 - Mouse shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	NM_010928 , NM_010928.1 , NM_010928.2 , BC059256 , BC172642
UniProt ID:	O35516
Summary:	Functions as a receptor for membrane-bound ligands Jagged-1 (JAG1), Jagged-2 (JAG2) and Delta-1 (DLL1) to regulate cell-fate determination (PubMed:10393120). Upon ligand activation through the released notch intracellular domain (NICD) it forms a transcriptional activator complex with RBPJ/RBPSUH and activates genes of the enhancer of split locus (PubMed:10393120, PubMed:18710934). Affects the implementation of differentiation, proliferation and apoptotic programs (PubMed:10393120, PubMed:18710934). May play an essential role in postimplantation development, probably in some aspect of cell specification and/or differentiation (By similarity). In collaboration with RELA/p65 enhances NFATc1 promoter activity and positively regulates RANKL-induced osteoclast differentiation (PubMed:18710934). Positively regulates self-renewal of liver cancer cells (By similarity). [UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).