

Product datasheet for **TL501616**

Pea15a Mouse shRNA Plasmid (Locus ID 18611)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Pea15a Mouse shRNA Plasmid (Locus ID 18611)
Locus ID:	18611
Synonyms:	Mat; Mat1; Pea; PEA-; PEA-15; Pea15; Pkcs; Pkcs15
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Pea15a - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 18611). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC038282 , NM_011063 , NM_011063.1 , NM_011063.2 , NM_011063.3 , BC092079
UniProt ID:	Q62048
Summary:	This gene encodes an adaptor protein that functions as a negative regulator of apoptosis induced by tumor necrosis factor-alpha, tumor necrosis factor-related apoptosis-inducing ligand, and Fas, through its interaction with fas-associated protein with death domain and caspase-8. It also regulates proliferation signaling by relocating the extracellular signal-regulated protein kinases 1 and 2 to the cytosol. The protein encoded by this gene contains an N-terminal death effector domain and a long, flexible C-terminal tail. In humans, the encoded protein is an endogenous substrate for protein kinase C. This protein is overexpressed in type 2 diabetes mellitus, where it may contribute to insulin resistance in glucose uptake. Multiple pseudogenes of this gene have been identified. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2016]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).