

Product datasheet for **TL501651V**

Pitx1 Mouse shRNA Lentiviral Particle (Locus ID 18740)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Pitx1 Mouse shRNA Lentiviral Particle (Locus ID 18740)
Locus ID:	18740
Synonyms:	Bft; P-OTX; Potx; Ptx1
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	Pitx1 - Mouse shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	BC012696 , NM_011097 , NM_011097.1 , NM_011097.2
UniProt ID:	P70314
Summary:	Sequence-specific transcription factor that binds gene promoters and activates their transcription. May play a role in the development of anterior structures, and in particular, the brain and facies and in specifying the identity or structure of hindlimb. Can independently activate and synergize with PIT-1 on pituitary-specific target gene promoters, thus may subserve functions in generating both precursor and specific cell phenotypes in the anterior pituitary gland and in several other organs. Can activate pituitary transcription of the proopiomelanocortin gene.[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).