

## Product datasheet for **TL501800V**

### Ptpa Mouse shRNA Lentiviral Particle (Locus ID 19262)

#### Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Ptpa Mouse shRNA Lentiviral Particle (Locus ID 19262)
Locus ID:	19262
Synonyms:	Ptpa; Ptpalpha; Rptalpha; Rptpa; Rptralpha
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	Ptpa - Mouse shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 <sup>7</sup> TU/ml.
RefSeq:	<a href="#">BC050930</a> , <a href="#">BC110683</a> , <a href="#">NM_001163688</a> , <a href="#">NM_008980</a> , <a href="#">NM_001355161</a> , <a href="#">NM_001163688.1</a> , <a href="#">NM_008980.1</a> , <a href="#">NM_008980.2</a>
UniProt ID:	<a href="#">P18052</a>
Summary:	Tyrosine protein phosphatase which is involved in integrin-mediated focal adhesion formation (PubMed:22801373). Following integrin engagement, specifically recruits BCAR3, BCAR1 and CRK to focal adhesions thereby promoting SRC-mediated phosphorylation of BRAC1 and the subsequent activation of PAK and small GTPase RAC1 and CDC42 (PubMed:22801373).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).