

## Product datasheet for **TL501825**

### Rab12 Mouse shRNA Plasmid (Locus ID 19328)

#### Product data:

|                           |   |
|---------------------------|---|
| Product Type:             | shRNA Plasmids  |
| Product Name:             | Rab12 Mouse shRNA Plasmid (Locus ID 19328)  |
| Locus ID:                 | 19328   |
| Synonyms:                 | 2900054P15Rik; AW555448; C77700   |
| Vector:                   | pGFP-C-shLenti (TR30023)  |
| E. coli Selection:        | Chloramphenicol (34 ug/ml)  |
| Mammalian Cell Selection: | Puromycin   |
| Format:                   | Lentiviral plasmids   |
| Components:               | Rab12 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 19328).<br>5µg purified plasmid DNA per construct<br>29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.  |
| RefSeq:                   | <a href="#">NM_024448</a> , <a href="#">NM_024448.1</a> , <a href="#">NM_024448.2</a> , <a href="#">BC038883</a> , <a href="#">BC065420</a> , <a href="#">BC152766</a>  |
| UniProt ID:               | <a href="#">P35283</a>  |
| Summary:                  | The small GTPases Rab are key regulators of intracellular membrane trafficking, from the formation of transport vesicles to their fusion with membranes. Rabs cycle between an inactive GDP-bound form and an active GTP-bound form that is able to recruit to membranes different set of downstream effectors directly responsible for vesicle formation, movement, tethering and fusion (By similarity). That Rab may play a role in protein transport from recycling endosomes to lysosomes regulating, for instance, the degradation of the transferrin receptor. Involved in autophagy.[UniProtKB/Swiss-Prot Function] |
| shRNA Design:             | These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .  |



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).