

## Product datasheet for **TL501866**

### Rbbp6 Mouse shRNA Plasmid (Locus ID 19647)

#### Product data:

|                           |   |
|---------------------------|---|
| Product Type:             | shRNA Plasmids  |
| Product Name:             | Rbbp6 Mouse shRNA Plasmid (Locus ID 19647)  |
| Locus ID:                 | 19647   |
| Synonyms:                 | 4933422O15Rik; AI316869; BB233631; C030034J04Rik; C77662; P2P-R; PACT   |
| Vector:                   | pGFP-C-shLenti (TR30023)  |
| E. coli Selection:        | Chloramphenicol (34 ug/ml)  |
| Mammalian Cell Selection: | Puromycin   |
| Format:                   | Lentiviral plasmids   |
| Components:               | Rbbp6 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 19647).<br>5µg purified plasmid DNA per construct<br>29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.  |
| RefSeq:                   | <a href="#">NM_011247</a> , <a href="#">NM_175023</a> , <a href="#">NM_011247.1</a> , <a href="#">NM_011247.2</a> , <a href="#">NM_175023.1</a> , <a href="#">NM_175023.2</a> ,<br><a href="#">NM_175023.3</a> , <a href="#">BC025874</a> , <a href="#">BC029229</a> , <a href="#">BC040426</a> , <a href="#">BC052528</a>  |
| UniProt ID:               | <a href="#">P97868</a>  |
| Summary:                  | E3 ubiquitin-protein ligase which promotes ubiquitination of YBX1, leading to its degradation by the proteasome (By similarity). May play a role as a scaffold protein to promote the assembly of the p53/TP53-MDM2 complex, resulting in increase of MDM2-mediated ubiquitination and degradation of p53/TP53; may function as negative regulator of p53/TP53, leading to both apoptosis and cell growth retardation (PubMed:17470788). Regulates DNA-replication and common fragile sites (CFS) stability in a ZBTB38- and MCM10-dependent manner. Controls ZBTB38 protein stability and abundance via ubiquitination and proteasomal degradation, and ZBTB38 in turn negatively regulates the expression of MCM10 which plays an important role in DNA-replication (PubMed:24726359).[UniProtKB/Swiss-Prot Function] |
| shRNA Design:             | These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .  |



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).