

Product datasheet for **TL501952**

Rtn2 Mouse shRNA Plasmid (Locus ID 20167)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Rtn2 Mouse shRNA Plasmid (Locus ID 20167)
Locus ID:	20167
Synonyms:	MMS10-P; Ms10p; Nspl1; NSPLI
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Rtn2 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 20167). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC031370 , NM_001025364 , NM_013648 , NM_013648.1 , NM_013648.2 , NM_013648.3 , NM_013648.4 , NM_013648.5 , NM_013648.6 , NM_001025364.1 , NM_001025364.2 , NM_001025364.3
UniProt ID:	O70622
Summary:	Inhibits amyloid precursor protein processing, probably by blocking BACE1 activity (By similarity). Enhances trafficking of the glutamate transporter SLC1A1/EAAC1 from the endoplasmic reticulum to the cell surface (By similarity). Plays a role in the translocation of SLC2A4/GLUT4 from intracellular membranes to the cell membrane which facilitates the uptake of glucose into the cell (PubMed:19720795).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).