

## Product datasheet for **TL501987V**

### Ccl2 Mouse shRNA Lentiviral Particle (Locus ID 20296)

#### Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Ccl2 Mouse shRNA Lentiviral Particle (Locus ID 20296)
Locus ID:	20296
Synonyms:	AI323594; HC11; JE; MCA; MCAF; MCP; MCP-; MCP-1; MCP1; Scy; Scya2; Sig; Sigje; SMC-C; SMC-CF
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	Ccl2 - Mouse shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 <sup>7</sup> TU/ml.
RefSeq:	<a href="#">NM_011333</a> , <a href="#">NM_011333.1</a> , <a href="#">NM_011333.2</a> , <a href="#">NM_011333.3</a> , <a href="#">BC145867</a> , <a href="#">BC055070</a> , <a href="#">BC145869</a>
UniProt ID:	<a href="#">P10148</a>
Summary:	This gene is one of several cytokine genes clustered on chromosome 11. Chemokines are a superfamily of secreted proteins involved in immunoregulatory and inflammatory processes. The superfamily is divided into four subfamilies based on the arrangement of N-terminal cysteine residues of the mature peptide. This chemokine is a member of the CC subfamily which is characterized by two adjacent cysteine residues. This cytokine displays chemotactic activity for monocytes and memory T cells but not for neutrophils. The human ortholog has been implicated in the pathogenesis of diseases characterized by monocytic infiltrates, such as psoriasis, rheumatoid arthritis, and atherosclerosis. [provided by RefSeq, Sep 2015]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).