

Product datasheet for **TL502041**

Six6 Mouse shRNA Plasmid (Locus ID 20476)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Six6 Mouse shRNA Plasmid (Locus ID 20476)
Locus ID:	20476
Synonyms:	Opt; Optx2; Si; Six9
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Six6 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 20476). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_011384 , NM_011384.1 , NM_011384.2 , NM_011384.3 , NM_011384.4 , NM_011384.5 , BC138838 , BC138839
UniProt ID:	Q9QZ28
Summary:	This gene encodes a homeobox protein that is similar to the Drosophila 'sine oculis' gene product. This gene is found in a cluster of related genes on chromosome 12 and is thought to be involved in eye development. The encoded transcription factor regulates early progenitor cell proliferation during mammalian retinogenesis and pituitary development. Mice lacking this gene exhibit abnormal development of the suprachiasmatic nucleus and circadian rhythms. [provided by RefSeq, Sep 2015]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).