

## Product datasheet for **TL502084V**

### Sumo3 Mouse shRNA Lentiviral Particle (Locus ID 20610)

#### Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Sumo3 Mouse shRNA Lentiviral Particle (Locus ID 20610)
Locus ID:	20610
Synonyms:	2810014B19Rik; D10Ertd345; D10Ertd345e; SMT; SMT3A; Smt3h; Smt3h1; SUMO-3
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	Sumo3 - Mouse shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 <sup>7</sup> TU/ml.
RefSeq:	<a href="#">BC115488</a> , <a href="#">BC115489</a> , <a href="#">NM_001301671</a> , <a href="#">NM_001301672</a> , <a href="#">NM_001301673</a> , <a href="#">NM_019929</a> , <a href="#">NR_125903</a> , <a href="#">NR_125904</a> , <a href="#">NM_019929.1</a> , <a href="#">NM_019929.2</a> , <a href="#">NM_019929.3</a> , <a href="#">NM_019929.4</a> , <a href="#">NM_001301672.1</a> , <a href="#">NM_001301671.1</a> , <a href="#">NM_001301673.1</a>
UniProt ID:	<a href="#">Q9Z172</a>
Summary:	This gene encodes a member of the small ubiquitin-like modifier family. The encoded protein may regulate a variety of proteins in many pathways via a post-translational modification, known as SUMOylation. This activity may play a role in a wide variety of cellular processes, including nuclear transport, DNA replication and repair, mitosis, transcriptional regulation, and signal transduction. Disruption of some of these processes has been associated with cerebral ischemia, neural dysfunction, and heart disease. A pseudogene of this gene has been defined on the X chromosome. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2014]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).