

Product datasheet for **TL502176V**

Suc1g2 Mouse shRNA Lentiviral Particle (Locus ID 20917)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Suc1g2 Mouse shRNA Lentiviral Particle (Locus ID 20917)
Locus ID:	20917
Synonyms:	AF171077; AW556404; D6Wsu120e; G-SCS; GTPSCS; SCS-betaG
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	Suc1g2 - Mouse shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	BC080781 , NM_001326558 , NM_011507 , NM_011507.1 , NM_011507.2 , NM_011507.3 , BC030947 , BC043312 , BC054425
UniProt ID:	Q9Z2I8
Summary:	GTP-specific succinyl-CoA synthetase functions in the citric acid cycle (TCA), coupling the hydrolysis of succinyl-CoA to the synthesis of GTP and thus represents the only step of substrate-level phosphorylation in the TCA. The beta subunit provides nucleotide specificity of the enzyme and binds the substrate succinate, while the binding sites for coenzyme A and phosphate are found in the alpha subunit.[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).