

Product datasheet for **TL502234V**

Tfap2b Mouse shRNA Lentiviral Particle (Locus ID 21419)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Tfap2b Mouse shRNA Lentiviral Particle (Locus ID 21419)
Locus ID:	21419
Synonyms:	AI606113; AP-2(beta); E130018K07Rik; Tcfap2b
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	Tcfap2b - Mouse shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	<u>NM_001025305</u> , <u>NM_001286340</u> , <u>NM_009334</u> , <u>NM_009334.1</u> , <u>NM_009334.2</u> , <u>NM_009334.3</u> , <u>NM_001025305.1</u> , <u>NM_001025305.2</u> , <u>NM_001286340.1</u> , <u>BC146412</u> , <u>BC148852</u>
UniProt ID:	<u>Q61313</u>
Summary:	Sequence-specific DNA-binding protein that interacts with inducible viral and cellular enhancer elements to regulate transcription of selected genes. AP-2 factors bind to the consensus sequence 5'-GCCNNNGGC-3' and activate genes involved in a large spectrum of important biological functions including proper eye, face, body wall, limb and neural tube development. They also suppress a number of genes including MCAM/MUC18, C/EBP alpha and MYC. AP-2-beta appears to be required for normal face and limb development and for proper terminal differentiation and function of renal tubular epithelia.[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).