

Product datasheet for **TL502237**

Tfec Mouse shRNA Plasmid (Locus ID 21426)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Tfec Mouse shRNA Plasmid (Locus ID 21426)
Locus ID:	21426
Synonyms:	BB107417; bHLHe34; Tcfec
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Tfec - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 21426). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC098494 , NM_031198 , NM_031198.1 , NM_031198.2 , NM_031198.3 , BC052832 , BC064770
UniProt ID:	Q9WTW4
Summary:	Transcriptional regulator that acts as a repressor or an activator. Acts as a transcriptional transactivator on the proximal promoter region of the tartrate-resistant acid phosphatase (TRAP) E-box containing promoter. Collaborates with MITF in target gene activation. Acts as a transcriptional repressor on minimal promoter containing element F (that includes an E-box sequence) (By similarity). Binds to element F in an E-box sequence-specific manner (By similarity). Acts as a transcriptional repressor on minimal promoter containing mu E3 enhancer sequence (By similarity). Binds to mu E3 DNA sequence of the immunoglobulin heavy-chain gene enhancer (By similarity). Binds DNA in a homo- or heterodimeric form. [UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).