

## **Product datasheet for TL502347V**

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

### **Tub Mouse shRNA Lentiviral Particle (Locus ID 22141)**

#### **Product data:**

**Product Type:** shRNA Lentiviral Particles

**Product Name:** Tub Mouse shRNA Lentiviral Particle (Locus ID 22141)

Locus ID: 22141 Synonyms: rd5

**Vector:** pGFP-C-shLenti (TR30023)

Format: Lentiviral particles

Components: Tub - Mouse shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble

control), 0.5 ml each, >10^7 TU/ml.

**RefSeq:** <u>BC096554</u>, <u>NM 021885</u>, <u>NM 021885.1</u>, <u>NM 021885.2</u>, <u>NM 021885.3</u>, <u>NM 021885.4</u>, <u>BC030921</u>

UniProt ID: P50586

**Summary:** Functions in signal transduction from heterotrimeric G protein-coupled receptors. Binds to

membranes containing phosphatidylinositol 4,5-bisphosphate. Can bind DNA (in vitro). May

contribute to the regulation of transcription in the nucleus. Could be involved in the hypothalamic regulation of body weight. Contribute to stimulation of phagocytosis of apoptotic retinal pigment epithelium (RPE) cells and macrophages (By similarity).

[UniProtKB/Swiss-Prot Function]

shRNA Design: These shRNA constructs were designed against multiple splice variants at this gene locus. To

be certain that your variant of interest is targeted, please contact <u>techsupport@origene.com</u>. If you need a special design or shRNA sequence, please utilize our custom shRNA service.





# Performance Guaranteed:

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).