

Product datasheet for **TL502377V**

Dpysl3 Mouse shRNA Lentiviral Particle (Locus ID 22240)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Dpysl3 Mouse shRNA Lentiviral Particle (Locus ID 22240)
Locus ID:	22240
Synonyms:	CRMP; CRMP-4; DRP-3; TUC4; U; Ul; Ulip; ULIP-1; Ulip1
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	Dpysl3 - Mouse shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	BC023003 , NM_001136086 , NM_001291455 , NM_009468 , NM_009468.1 , NM_009468.2 , NM_009468.3 , NM_009468.4 , NM_009468.5 , NM_001136086.1 , NM_001136086.2 , NM_001291455.1
UniProt ID:	Q62188
Summary:	This gene encodes a protein that belongs to the TUC (TOAD-64/Ulip/CRMP) family of proteins. Members of this family are phosphoproteins that function in axonal guidance and neuronal differentiation during development and regeneration of the nervous system. A mutation in the human gene is associated with amyotrophic lateral sclerosis. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Apr 2014]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).